

☎ 86-20-83565714

☎ 86-13005107206

✉ info@lingchendental.com

📘 www.facebook.com/lingchendental

Web www.lingchendental.com

📍 Room3304,Jinying Building, No.316 Huan Shi
Zhong Road,Yuexiu Area,Guangzhou,China

ALL RIGHTS RESERVED

LINGCHEN

LINGCHEN DENTAL

Dental Chair / Autoclave / Simulator

Manufacturer &
One-stop Dental Solution



Premium
Quality




Reliable
Performance



Global
Trust



Complete
Solution



©2026 LINGCHEN

COMPANY FILE

Dental Chair / Autoclave / Simulator Manufacturer & One-stop Dental Solution

-Founded in 2009, located in Guangzhou China.

Specialized in designing, developing, producing and selling dental equipment, targeting the Middle East, Africa, and Latin America with 16 years experience.


-Operating 3 production lines with 40+ employees in a 2000+ square meter factory.

-Export to 70+ countries, supporting 30+ universities.

Certificates



Project



Exhibition



LINGCHEN

Implant Dental Chair TAOS1800i

Simplify Implant Procedures, Eliminate Headaches

- Save workspace, no need additional trolley.
- Helps save time before and after surgery more than 30mins.
- 3 trays allow to put the tools in good regulation.

Built-in
Implant Motor

The First
Real Implant
Dental Chair



Tray 2
For Implant Accessories

Tray 1
For Surgery Tools

Tray 3
For Implant Motor Things



King Dental Chair TAOS1800



Full options

LONGCHEN



Luna high speed handpiece 2pcs

Full set low speed handpiece

Lingchen built-in LED scaler

Q7 curing light

CCD oral camera with monitor

Air compressor AC1100W



Touch screen + button control

New distiller water system

Sensor water giving

8 eyes LED lamp

Customer Feedback



Prince Dental Chair TAOS800L


- Luxury injection molding operation tray.
- 2 meters long cushion.
- Vacuum cleaning system.

Option: electric motor



Full options

LONGCHEN



Dental Center Clinic Unit *TAOS1800c*

Used in Surgery, RCT, Implant, Education.

- A imaginative design - saving time and cost.
- Providing a large operating space, which is needed for all the operations and treatment.




Queen Dental Chair *TAOS900*

LINGCHEN

Built-in camera



Fortiva Plus Dental Chair LC900



Details



Full options LINGCHEN

	
High Speed Handpiece	Low Speed Handpiece
	
Curing Light	Built-in Scaler
	
CCD Intra Oral Camera	Air Compressor AC1100

Q1 Kids Dental Chair

-The toy car design makes it easy to build a children's trust and reduce their fear.



Q2 Kids Dental Chair

-Tom and Jerry International language for kids to get their interest.



Air Compressor & Vacuum Pump



AC1100
1100W, 50L



AC1500 with air dryer
1500W, 50L



AC1500 with cabinet



Vacuum Pump 750W
Used for two dental chair

LONGCHEN



AC1500-3 with air dryer
3300W, 150L




AC1500-5 with air dryer
5500W, 200L

TS18 Dental Autoclave



- Class B
- 22 mins whole circle
- Accessories after-sales support

Sterilization Cycle



- 22 mins to finish autoclave circle
- Euro Technology -Class B, apply for metal, fabric, cotton, plastic (internal & external)
- Digital LED display, 121°C/134°C

LINGCHEN



Project



Spare Part



Factory



Feedback

Athena Dental Simulator SS01

Usage:

- Save space & cost, 30 square meter can put 15 units.
- Operation tray design same as dental chair,help student adopt the dental chair using by fast way.



Options 1



Camera system which can show and make video from teacher, improving the efficiency between teaching and practicing.

Chip: Sony 4K@30fps, 1080P,
30 frames 8 million pixels
Effective pixels: 3840(H) x 2160(V)
Focus mode: auto focus (20x zoom)

Options 2



Goggle loupe
3.5X magnificent
with LED light

Options 3



Compressor 550W
22L tank



LONGCHEN

Luban Dental Simulator SS03

- Small separated package, to save the shipping cost and easy to move into the room.
- Observation system two ways, improving the efficiency between teaching and practising.
- Cabinet can be customized.
- Designed as the lab bench, 1 table 2 usage and the table surface is fireproof, more safety.



Jaw Model
(Order No.:1154)



Pedo Jaw Model
(Order No.:1155)



Adult Teeth
(Order No.:1156)



Kids Teeth
(Order No.:1156-1)



Option 1: observation system

Teacher+ Student HD video terminal:

Main Server + Teacher management system
+ Control Pad + Teacher HD camera
+ Student HD camera



Teacher

Photos & Videos



Students

Which can share the photos and videos between the teachers and students

Option 2: microscope



- The continuous adjustable range of varifocus objective is 300mm.
- Continue zooming (Manual)
- 3.36X~20X (with 10X Eyepiece)
- Stepless zoom with seamless clarity across the entire range
- **Adjustable range:** 55 -75mm
- Integrated 4K imaging module (4K/1080P adjustable)
- Supports video recording and photo capture
- Supports USB drive storage

LONGCHEN

Apollo Dental Simulator SS02

-Simulation system can be used in multimedia education and false tooth repair, which is necessary in stomatology college.

1 unit 2 students

- Both sides are available for students to practice.
- Save space and shipping costs

- Phantom Head*2
- Hanging Turbine*2
- LED Lamp*2
- Dentist Stool*2
- Drawer*2



Structural Diagram

